



Essential Nano LED Light

Technical Specs

Control Unit	Focus™ Xtend
Lamp weight	9 g
Lamp material	Aluminium
Operating time	6-24 hours depending on chosen light intensity
Charging time	Approx. 2.5 hours
LED diode	Natural white (6,500 Kelvin)
Maximum light intensity at 25 cm	40,000 lux
Light intensity at working distance 35 cm	18,000 lux
Light spot at working distance 35 cm	80 mm
Brightness setting	High, Medium and Low
Connection system	Focus™ Connect System
Cable type	Grey super flex
Warranty	3 years on Control Unit Xtend, Battery Xtend, lamp and accessories. 1 year on detachable cable and backplate.