



Essential One LED Light

Technical Specs

Control Unit	Essential One, Wireless
Lamp weight	7 g (battery weight is 16 g)
Lamp material	Aluminium
Operating time	Up to 2 hours per battery (3 batteries included)
Charging time	Approx. 1.5 hours
LED diode	Natural white (5,750 Kelvin)
Maximum light intensity at 25 cm	33,000 lux
Light intensity at working distance 35 cm	14,000 lux
Light spot at working distance 35 cm	65 mm
Brightness setting	2: High and Low
Connection system	Essential One framepart
Warranty	Battery 1 year/Lamp 2 years